



# Precision Angled Screw Restoration Overview



## PRECISION ANGLED SCREW RESTORATION OVERVIEW

THERE ARE FOUR STEPS TO CREATING AND DELIVERING A PRECISION ANGLED SCREW RESTORATION. THE STEPS ARE AS FOLLOWS:


### 1 DIGITIZE THE IMPRESSION

There are two primary ways to create a digital impression.

**Option A:** The first and most common method is to take an implant level impression, pour a stone model, place a software-dependent scan body onto the implant analog and scan the model using a 3D tabletop digital scanner.

For implant level impressions, refer to the **open tray technique using the direct pick-up coping** module, the **closed tray technique using the indirect transfer coping** module and the **closed tray pick-up technique using the snap coping** module.

**Option B:** The second method is to take an intra-oral digital impression by placing a software-dependent scan body into the implant and scan the scan body and surrounding dentition using a handheld 3D scanner.



**Important:**

In order to design a Precision Angled Screw restoration, a BioHorizons scan body must be used and the BioHorizons CAD/CAM Abutments & Model library must be downloaded from the Vulcan website. <http://www.vulcandental.com/Home/Downloads>


### 2 DESIGN THE CROWN

The file that is created during the digital impression is imported into the design software that will be used by the technician to design the angled screw channel restoration.



**Note:**

The precision angled screw and driver are compatible with the BioHorizons Hybrid Base abutments and Tissue-Level Hybrid Base abutments.



## PRECISION ANGLED SCREW RESTORATION OVERVIEW

CREATING AND DELIVERING A PRECISION ANGLED SCREW CHANNEL RESTORATION ON A BIOHORIZONS HYBRID BASE ABUTMENT. STEPS CONTINUED:

### 3 CEMENT THE CROWN

The laboratory will use the Hybrid Base abutment and stone or printed model to complete the crown following normal laboratory procedures. The gold Precision Angled Screw should remain in the Hybrid Base abutment during the cementation of the crown to the abutment.



**Important:**

During cementation of the crown to the abutment, prevent cement from entering into the angled screw channel of the crown and the internal aspect of the abutment. The final angled screw channel restoration will have a captured screw within the assembly.



### 4 DELIVER THE FINAL RESTORATION

The final restoration should be sanitized following standard clinical procedures. When utilizing the angled screw channel feature, **ONLY** a Precision Angled Driver should be used to tighten the gold Precision Angled Screw that is captured in the restoration.

Secure the restoration onto the implant using the Precision Angled Driver. Hand tighten.

Tighten the abutment screw to 30Ncm using the Precision Angled Driver and a calibrated torque wrench.



**Important:**

When using the calibrated torque wrench, apply a steady, downward pressure along the axis of the driver to maintain optimum driver and screw engagement.





2300 Riverchase Center, Suite 825  
Birmingham, AL 35244  
Toll free: 1.844.484.2301 • Direct: 205.484.2301  
[www.vulcandental.com](http://www.vulcandental.com)