

Vulcan Library Supported Elos Accurate® Products		IO Scan Body	Model Analogs	Custom Titanium	Hybrid Base Engaging		Hybrid Base Non-Engaging		Prosthetic screw for straight channels	Prosthetic screw - hexalobular	Torque	
Implant System	Platform	IO	PMA	TIABUT	HBE	HBE-H	HBN	HBN-H	AS	Screwdriver:	AS hexalobular	Ncm:
BioHorizons Internal	3.0	IO 11A-A		TIABUT101								30 Ncm
	3.5	IO 11A-B		TIABUT102								
	4.5	IO 11A-C		TIABUT103								
	5.7	IO 11A-D		TIABUT104								
BioHorizons Multi-Unit Abutment	MUA	IO 11B-A										15 Ncm
Camlog® Bar Abutment	3.3-4.3	IO 9B-A	PMA-CBA43-1				HBN-CBA43-1				AS-CBAM1604A-1	15 Ncm
Astra Tech Implant System TX	3.0	IO 3A-A	PMA-ATO30-1		HBE-ATO30-1				AS-ATOM148POS	"0.050" Hex screwdriver	AS-ATOM1408A-1	15 Ncm
	3.5/4.0	IO 3A-B	PMA-ATO35-1	TIABUT315	HBE-ATO35/40-1	HBE-ATO35/40H-1			AS-ATOM168POS		AS-ATOM1608A-1	20 Ncm
	4.5/5.0	IO 3A-C	PMA-ATO45-1	TIABUT316	HBE-ATO45/50-1	HBE-ATO45/50H-1			AS-ATOM2010POS		AS-ATOM2010A-1	25 Ncm
Astra Tech Implant System EV	3.0	IO 3B-A	PMA-ATE30-1		HBE-ATE30-1				AS-ATEM148POS	"0.050" Hex screwdriver	AS-ATEM1408A-1	25 Ncm
	3.6	IO 3B-B	PMA-ATE36-1	TIABUT311	HBE-ATE36-1	HBE-ATE36H-1			AS-ATEM168POS		AS-ATEM1608A-1	
	4.2	IO 3B-C	PMA-ATE42-1	TIABUT312	HBE-ATE42-1	HBE-ATE42H-1			AS-ATEM188POS		AS-ATEM1808A-1	
	4.8	IO 3B-D	PMA-ATE48-1	TIABUT313	HBE-ATE48-1	HBE-ATE48H-1			AS-ATEM208POS		AS-ATEM2008A-1	
Astra Tech Implant System UniAbutment EV	UA33°	IO 3D-A	PMA-AUE33-1				HBN-AUE33-1	HBN-AUE33H-1			AS-AUEM1805A-1	15 Ncm
Brånemark System	3.3/3.5 NP	IO 6A-A	PMA-BRA35-1	TIABUT317					AS-BRAM167POS	Nobel Biocare Unigrip		35 Ncm
	3.75/4.0 RP	IO 6A-B	PMA-BRA41-1	TIABUT318					AS-BRAM207POS			
	5.0/5.5 WP	IO 6A-C	PMA-BRA51-1	TIABUT319					AS-BRAM257POS			
Nobel Biocare® Conical Connection and NobelActive®	3.0	IO 2B-D / SA	PMA-NBA30-1	TIABUT303	HBE-NBA30-1				AS-NBAM1407POS	Nobel Biocare Unigrip	AS-NBAM1407A-1	15 Ncm
	3.5/3.75 NP	IO 2B-A / SA	PMA-NBA35-1	TIABUT304	HBE-NBA35-1	HBE-NBA35H-1	HBN-NBA35-1	HBN-NBA35H-1	AS-NBAM168POS		AS-NBAM1608A-1	35 Ncm
	4.3/5.0 RP	IO 2B-B / SA	PMA-NBA43-1	TIABUT305	HBE-NBA43-1	HBE-NBA43H-1	HBN-NBA43/50-1	HBN-NBA43/50H-1	AS-NBAM207POS		AS-NBAM2007A-1	
	5.5 WP	IO 2B-C / SA	PMA-NBA60-1	TIABUT306	HBE-NBA60-1	HBE-NBA60H-1	HBN-NBA60-1					
Nobel Biocare® Multi-Unit Abutment	NP/RP	IO 2C-A	PMA-MUA45-1				HBN-MUA45-1	HBN-MUA45H-1	AS-MUAM1403POS	Nobel Biocare Unigrip	AS-MUAM1403A-1	15 Ncm
	WP	IO 2C-B	PMA-MUA60-1				HBN-MUA60-1		AS-MUAM1804POS		AS-MUAM1804A-1	
NobelReplace®	3.5 NP	IO 2A-A	PMA-NBR35-1	TIABUT307	HBE-NBR35-1	HBE-NBR35H-1	HBN-NBR35-1	HBN-NBR35H-1	AS-NBRM1808POS	Nobel Biocare Unigrip	AS-NBRM1808A-1	35 Ncm
	4.3 RP	IO 2A-B	PMA-NBR43-1	TIABUT308	HBE-NBR43-1	HBE-NBR43H-1	HBN-NBR43-1	HBN-NBR43H-1	AS-NBRM2010POS		AS-NBRM2010A-1	
	5.0 WP	IO 2A-C	PMA-NBR50-1	TIABUT309	HBE-NBR50-1		HBN-NBR50-1					
	6.0	IO 2A-D	PMA-NBR60-1	TIABUT310	HBE-NBR60-1		HBN-NBR60-1					
Straumann® Bone Level	3.3 NC	IO 4A-A	PMA-SBO33-1	TIABUT301	HBE-SBO33-1	HBE-SBO33H-1	HBN-SBO33-1		AS-SBOM1608POS	SCS screwdriver	AS-SBOM1608A-1	35 Ncm
	4.1/4.8 RC	IO 4A-B	PMA-SBO41-1	TIABUT302	HBE-SBO41/48-1	HBE-SBO41/48H-1	HBN-SBO41/48-1					
Straumann® Screw-retained Abutment	3.5	IO 4D-A	PMA-SSA35-1				HBN-SSA35-1	HBN-SSA35H-1			AS-SSAM1404A-1	15 Ncm
	4.6	IO 4D-B	PMA-SSA46-1				HBN-SSA46-1	HBN-SSA46H-1				
Price		\$74	\$37	\$169	\$74	\$74	\$74	\$74	\$31		\$37	