

Elos Accurate®		IO Scan Body	Model Analogs	Custom Titanium	Hybrid Base Engaging		Hybrid Base Non-Engaging		Prosthetic screw for straight channels		Prosthetic screw - hexalobular	Torque	
Platform		IO	PMA	TIABUT	HBE	HBE-H	HBN	HBN-H	AS	Screwdriver	AS hexalobular	Ncm:	
BIOHORIZONS	BioHorizons Internal	3.0	IO 11A-A	PMA-BHI30-1	TIABUT101							30 Ncm	
		3.5	IO 11A-B	PMA-BHI35-1	TIABUT102								
		4.5	IO 11A-C	PMA-BHI45-1	TIABUT103								
		5.7	IO 11A-D	PMA-BHI57-1	TIABUT104								
	BioHorizons Multi-Unit Abutment	MUA	IO 11B-A	PMA-BHM45-1									
CAMLOG	Camlog®	3.3	IO 9A-A	PMA-CAS33-1	TIABUT201							20 Ncm	
		3.8	IO 9A-B	PMA-CAS38-1	TIABUT202								
		4.3	IO 9A-C	PMA-CAS43-1	TIABUT204								
		5,0	IO9A-D	PMA-CAS50-1	TIABUT206								
	6,0	IO9A-D	PMA-CAS60-1	TIABUT208									
	Camlog® Bar Abutment	3.3-4.3	IO 9B-A	PMA-CBA43-1			HBN-CBA43-1				AS-CBAM1604A-1	15 Ncm	
DENTSPLY SIRONA	Ankylos®	3.5-7.0	IO 8B-A	PMA-DAN35/70-1									
	Ankylos Balance Base Abutment	BB NP	IO 8C-A	PMA-DBA40-1			HBN-DBA40-1	HBN-DBA40H-1	AS-DBAM1605POS		AS-DBAM1605A-1	10 Ncm	
	Astra Tech Implant System TX	3.0	IO 3A-A	PMA-ATO30-1		HBE-ATO30-1				AS-ATOM148POS	0.050 Hex screwdriver	AS-ATOM1408A-1	15 Ncm
		3.5/4.0	IO 3A-B	PMA-ATO35-1	TIABUT315	HBE-ATO35/40-1	HBE-ATO35/40H-1			AS-ATOM168POS		AS-ATOM1608A-1	20 Ncm
		4.5/5.0	IO 3A-C	PMA-ATO45-1	TIABUT316	HBE-ATO45/50-1	HBE-ATO45/50H-1			AS-ATOM2010POS		AS-ATOM2010A-1	25 Ncm
	Astra Tech UniAbutment	20°	IO 3D-B	PMA-ATU20-1				HBN-ATU20-1				AS-ATUM1405A-1	15 Ncm
		45°	IO 3D-C	PMA-ATU45-1				HBN-ATU45-1					
	Astra Tech Implant System EV	3.0	IO 3B-A	PMA-ATE30-1		HBE-ATE30-1				AS-ATEM148POS	0.050 Hex screwdriver	AS-ATEM1408A-1	25 Ncm
		3.6	IO 3B-B	PMA-ATE36-1	TIABUT311	HBE-ATE36-1	HBE-ATE36H-1			AS-ATEM168POS		AS-ATEM1608A-1	
		4.2	IO 3B-C	PMA-ATE42-1	TIABUT312	HBE-ATE42-1	HBE-ATE42H-1			AS-ATEM188POS		AS-ATEM1808A-1	
		4.8	IO 3B-D	PMA-ATE48-1	TIABUT313	HBE-ATE48-1	HBE-ATE48H-1						
		5.4	IO 3B-E	PMA-ATE54-1	TIABUT314	HBE-ATE54-1				AS-ATEM208POS		AS-ATEM2008A-1	
Astra Tech Implant System Profile EV	4.2	IO 3C-A	PMA-ATE42P-1										
	4.8	IO 3C-B	PMA-ATE48P-1										
Astra Tech Implant System UniAbutment EV	UA33°	IO 3D-A	PMA-AUE33-1				HBN-AUE33-1	HBN-AUE33H-1			AS-AUEM1805A-1	15 Ncm	
Xive®	3.0	IO 8A-A	PMA-DXI30-1										
	3.4	IO 8A-B	PMA-DXI34-1										
	3.8	IO 8A-C	PMA-DXI38-1										
	4.5	IO 8A-D	PMA-DXI45-1										
	5,5	IO 8A-D	PMA-DXI55-1										
HIOSSEN	Hiossen ET	Mini - 3.2	IO 12A-A	PMA-HET32-1	TIABUT330				AS-HETM1610POS			20 Ncm	
		Mini - 3.5		PMA-HET35-1									
		Regular	IO 12A-B	PMA-HET45-1	TIABUT331				AS-HETM2008POS			30 Ncm	
NEOSS	Neoss® Access Abutment	5,5	IO 1B-A	PMA-NAC40-1									
	Neoss Implant System	3.25	IO 1A-B	PMA-NIS33-1									
		3.5-5.5	IO 1A-A	PMA-NIS35/55-1									
NOBEL BIOCARE	Brånemark System	3.3/3.5 NP	IO 6A-A	PMA-BRA35-1	TIABUT317	HBE-BRA35-1		HBN-BRA35-1		AS-BRAM167POS	Nobel Biocare Unigrip	AS-BRAM1607A-1	35 Ncm
		3.75/4.0 RP	IO 6A-B	PMA-BRA41-1	TIABUT318	HBE-BRA41-1	HBE-BRA41H-1	HBN-BRA41-1	HBN-BRA41H-1	AS-BRAM207POS		AS-BRAM2007A-1	
		5.0/5.5 WP	IO 6A-C	PMA-BRA51-1	TIABUT319	HBE-BRA51-1		HBN-BRA51-1		AS-BRAM257POS		AS-BRAM2507A-1	
	Nobel Biocare® Conical Connexion and NobelActive®	3.0	IO 2B-D	PMA-NBA30-1	TIABUT303	HBE-NBA30-1				AS-NBAM1407POS	Nobel Biocare Unigrip	AS-NBAM1407A-1	15 Ncm
		3.5/3.75 NP	IO 2B-A	PMA-NBA35-1	TIABUT304	HBE-NBA35-1	HBE-NBA35H-1	HBN-NBA35-1	HBN-NBA35H-1	AS-NBAM168POS		AS-NBAM1608A-1	
		4.3 RP	IO 2B-B	PMA-NBA43-1	TIABUT305	HBE-NBA43-1	HBE-NBA43H-1	HBN-NBA43/50-1	HBN-NBA43/50H-1	AS-NBAM207POS		AS-NBAM2007A-1	
	Nobel Biocare® Multi-Unit Abutment	NP/RP	IO 2C-A	PMA-MUA45-1				HBN-MUA45-1	HBN-MUA45H-1	AS-MUAM1403POS	Nobel Biocare Unigrip	AS-MUAM1403A-1	15 Ncm
		WP	IO 2C-B	PMA-MUA60-1				HBN-MUA60-1		AS-MUAM1804POS		AS-MUAM1804A-1	
	NobelReplace®	3.5 NP	IO 2A-A	PMA-NBR35-1	TIABUT307	HBE-NBR35-1	HBE-NBR35H-1	HBN-NBR35-1	HBN-NBR35H-1	AS-NBRM1808POS	Nobel Biocare Unigrip	AS-NBRM1808A-1	35 Ncm
		4.3 RP	IO 2A-B	PMA-NBR43-1	TIABUT308	HBE-NBR43-1	HBE-NBR43H-1	HBN-NBR43-1	HBN-NBR43H-1				
5.0 WP		IO 2A-C	PMA-NBR50-1	TIABUT309	HBE-NBR50-1		HBN-NBR50-1		AS-NBRM2010POS	AS-NBRM2010A-1			
6.0		IO 2A-D	PMA-NBR60-1	TIABUT310	HBE-NBR60-1		HBN-NBR60-1						
Price			\$74	\$39	\$179	\$74	\$74	\$74	\$74	\$19	\$25		

