

Elos Accurate®		Platform	IO Scan Body	Model Analogs	Model Analog Tool	Hybrid Base Engaging	Hybrid Base Non-Engaging	Prosthetic screw for straight channels	Prosthetic screw - hexalobular	Torque				
			IO	PMA	MA	AT	HBE	HBE-H	HBN	HBN-H	AS	Screwdriver:	AS hexalobular	Ncm:
BIOHORIZONS	BioHorizons Internal	3.0	IO 11A-A	PMA-BHI30-1										
		3.5	IO 11A-B	PMA-BHI35-1										
		4.5	IO 11A-C	PMA-BHI45-1										
		5.7	IO 11A-D	PMA-BHI57-1										
	BioHorizons Multi-Unit Abutment	MUA	IO 11B-A	PMA-BHM45-1										
CAMLOG	Camlog®	3.3	IO 9A-A	PMA-CAS33-1										
		3.8	IO 9A-B	PMA-CAS38-1	MA-CAS38-1	AT-16-1								
		4.3	IO 9A-C	PMA-CAS43-1	MA-CAS43-1									
		5.0	IO9A-D	PMA-CAS50-1	MA-CAS50-1									
	6.0	IO9A-D	PMA-CAS60-1											
	Camlog® Bar Abutment	3.3-4.3	IO 9B-A	PMA-CBA43-1	MA-CBA43-1	AT-16-1		HBN-CBA43-1				AS-CBAM1604A-1	15 Ncm	
DENTSPLY SIRONA	Ankylos®	3.5-7.0	IO 8B-A	PMA-DAN35/70-1	MA-DAN35/70-1	AT-18-1								
	Ankylos Balance Base Abutment	BB NP	IO 8C-A	PMA-DBA40-1	MA-DBA40-1	AT-16-1		HBN-DBA40-1	HBN-DBA40H-1	AS-DBAM1605POS		AS-DBAM1605A-1	10 Ncm	
	Astra Tech Implant System TX	3.0	IO 3A-A	PMA-ATO30-1	MA-ATO30-1	AT-14-1	HBE-ATO30-1				AS-ATOM148POS	"0.050" Hex screwdriver"	AS-ATOM1408A-1	15 Ncm
		3.5/4.0	IO 3A-B	PMA-ATO35-1	MA-ATO35-1	AT-16-1	HBE-ATO35/40-1	HBE-ATO35/40H-1			AS-ATOM168POS		AS-ATOM1608A-1	20 Ncm
		4.5/5.0	IO 3A-C	PMA-ATO45-1	MA-ATO45-1	AT-20-1	HBE-ATO45/50-1	HBE-ATO45/50H-1			AS-ATOM2010POS		AS-ATOM2010A-1	25 Ncm
	Astra Tech UniAbutment	20°	IO 3D-B	PMA-ATU20-1	MA-ATU20-1	AT-14-1			HBN-ATU20-1				AS-ATUM1405A-1	15 Ncm
		45°	IO 3D-C	PMA-ATU45-1	MA-ATU45-1				HBN-ATU45-1					
	Astra Tech Implant System EV	3.0	IO 3B-A	PMA-ATE30-1	MA-ATE30-1	AT-14-1	HBE-ATE30-1				AS-ATEM148POS	"0.050" Hex screwdriver"	AS-ATEM1408A-1	25 Ncm
		3.6	IO 3B-B	PMA-ATE36-1	MA-ATE36-1	AT-16-1	HBE-ATE36-1	HBE-ATE36H-1			AS-ATEM168POS		AS-ATEM1608A-1	
		4.2	IO 3B-C	PMA-ATE42-1	MA-ATE42-1	AT-18-1	HBE-ATE42-1	HBE-ATE42H-1			AS-ATEM188POS		AS-ATEM1808A-1	
4.8		IO 3B-D	PMA-ATE48-1	MA-ATE48-1	AT-20-1	HBE-ATE48-1	HBE-ATE48H-1			AS-ATEM208POS	AS-ATEM2008A-1			
5.4		IO 3B-E	PMA-ATE54-1	MA-ATE54-1			HBE-ATE54-1							
Astra Tech Implant System Profile EV	4.2	IO 3C-A	PMA-ATE42P-1	MA-ATE42P-1	AT-18-1									
	4.8	IO 3C-B	PMA-ATE48P-1	MA-ATE48P-1	AT-20-1									
Astra Tech Implant System UniAbutment EV	UA33°	IO 3D-A	PMA-AUE33-1	MA-AUE33-1	AT-18-1			HBN-AUE33-1	HBN-AUE33H-1			AS-AUEM1805A-1	15 Ncm	
Xive®	3.0	IO 8A-A	PMA-DXI30-1	MA-DXI30-1	AT-14-1									
	3.4	IO 8A-B	PMA-DXI34-1	MA-DXI34-1	AT-16-1									
	3.8	IO 8A-C	PMA-DXI38-1	MA-DXI38-1										
	4.5	IO 8A-D	PMA-DXI45-1	MA-DXI45-1										
	5.5	IO 8A-D	PMA-DXI55-1											
NEOSS	Neoss® Access Abutment		IO 1B-A	PMA-NAC40-1	MA-NAC40-1	AT-20-1								
	Neoss Implant System	3.25	IO 1A-B	PMA-NIS33-1	MA-NIS33-1	AT-16-1								
		3.5-5.5	IO 1A-A	PMA-NIS35/55-1	MA-NIS35/55-1	AT-20-1								
NOBEL BIO CARE	Brånemark System	3.3/3.5 NP	IO 6A-A	PMA-BRA35-1	MA-BRA35-1	AT-16-1	HBE-BRA35-1		HBN-BRA35-1		AS-BRAM167POS	Nobel Biocare Unigrip	AS-BRAM1607A-1	35 Ncm
		3.75/4.0 RP	IO 6A-B	PMA-BRA41-1	MA-BRA41-1	AT-20-1	HBE-BRA41-1	HBE-BRA41H-1	HBN-BRA41-1	HBN-BRA41H-1	AS-BRAM207POS		AS-BRAM2007A-1	
		5.0/5.5 WP	IO 6A-C	PMA-BRA51-1	MA-BRA51-1	AT-25-1	HBE-BRA51-1		HBN-BRA51-1		AS-BRAM257POS		AS-BRAM2507A-1	
	Nobel Biocare® Conical Connexion and NobelActive®	3.0	IO 2B-D / SA	PMA-NBA30-1	MA-NBA30-1	AT-14-1	HBE-NBA30-1				AS-NBAM1407POS	Nobel Biocare Unigrip	AS-NBAM1407A-1	15 Ncm
		3.5/3.75 NP	IO 2B-A / SA	PMA-NBA35-1	MA-NBA35-1	AT-16-1	HBE-NBA35-1	HBE-NBA35H-1	HBN-NBA35-1	HBN-NBA35H-1	AS-NBAM168POS		AS-NBAM1608A-1	
		4.3 RP	IO 2B-B / SA	PMA-NBA43-1	MA-NBA43-1	AT-20-1	HBE-NBA43-1	HBE-NBA43H-1	HBN-NBA43/50-1	HBN-NBA43/50H-1	AS-NBAM207POS		AS-NBAM2007A-1	
		5.0 RP			MA-NBA50-1									
	5.5 WP	IO 2B-C / SA	PMA-NBA60-1	MA-NBA60-1		HBE-NBA60-1	HBE-NBA60H-1	HBN-NBA60-1						
	Nobel Biocare® Multi-Unit Abutment	NP/RP	IO 2C-A	PMA-MUA45-1	MA-MUA45-1	AT-14-1			HBN-MUA45-1	HBN-MUA45H-1	AS-MUAM1403POS	Nobel Biocare Unigrip	AS-MUAM1403A-1	15 Ncm
		WP	IO 2C-B	PMA-MUA60-1	MA-MUA60-1	AT-18-1			HBN-MUA60-1		AS-MUAM1804POS		AS-MUAM1804A-1	
NobelReplace®	3.5 NP	IO 2A-A	PMA-NBR35-1	MA-NBR35-1	AT-18-1	HBE-NBR35-1	HBE-NBR35H-1	HBN-NBR35-1	HBN-NBR35H-1	AS-NBRM1808POS	Nobel Biocare Unigrip	AS-NBRM1808A-1	35 Ncm	
	4.3 RP	IO 2A-B	PMA-NBR43-1	MA-NBR43-1	AT-20-1	HBE-NBR43-1	HBE-NBR43H-1	HBN-NBR43-1	HBN-NBR43H-1	AS-NBRM2010POS				
	5.0 WP	IO 2A-C	PMA-NBR50-1	MA-NBR50-1			HBE-NBR50-1		HBN-NBR50-1					
	6.0	IO 2A-D	PMA-NBR60-1	MA-NBR60-1			HBE-NBR60-1		HBN-NBR60-1					

Price

\$74

\$37

\$37

\$32

\$74

\$74

\$74

\$74

\$31

\$37

Implant System	Platform	IO Scan Body	Model Analogs	Model Analog Tool	Hybrid Base Engaging		Hybrid Base Non-Engaging		Prosthetic screw for straight channels		Prosthetic screw - hexalobular	Torque	
		IO	PMA	MA	AT	HBE	HBE-H	HBN	HBN-H	AS	Screwdriver:	AS hexalobular	Ncm:
STRAUMANN	Straumann® Bone Level	3.3 NC	IO 4A-A	PMA-SBO33-1	AT-16-1	HBE-SBO33-1	HBE-SBO33H-1	HBN-SBO33-1		AS-SBOM1608POS	SCS screwdriver	AS-SBOM1608A-1	35 Ncm
		4.1 RC	IO 4A-B	PMA-SBO41-1		HBE-SBO41/48-1	HBE-SBO41/48H-1	HBN-SBO41/48-1					
		4.8 RC		MA-SBO48-1									
	Straumann® Screw-retained Abutment	3.5	IO 4D-A	PMA-SSA35-1				HBN-SSA35-1	HBN-SSA35H-1		AS-SSAM1404A-1	15 Ncm	
		4.6	IO 4D-B	PMA-SSA46-1				HBN-SSA46-1	HBN-SSA46H-1				
	Straumann® Standard and Standard Plus	3.5 NNC	IO 4C-A										
4.8 RN		IO 4B-A	PMA-SSY48-1	MA-SSY48-1	AT-20-1							35 Ncm	
6.5 WN		IO 4B-B	PMA-SSY65-1	MA-SSY65-1									
ZIMMER BIOMET	Certain® Internal Connection	3.4	IO 7B-A	PMA-BCE34-1		MA-BCE34-1	AT-16-1						
		4.1	IO 7B-B	PMA-BCE41-1	MA-BCE41-1								
		5.0		PMA-BCE50-1	MA-BCE50-1								
		6.0		PMA-BCE60-1	MA-BCE60-1								
External Hex Connection	3.4	IO 7A-A	PMA-BEH34-1	MA-BEH34-1	AT-20-1								
	4.1	IO 7A-B	PMA-BEH41-1	MA-BEH41-1									
	5.0		PMA-BEH50-1	MA-BEH50-1									
	6.0		PMA-BEH60-1	MA-BEH60-1									
Tapered Screw-Vent® Implant	3.5	IO 5A-A	PMA-ZSV35-1	MA-ZSV35-1	AT-18-1						30 Ncm		
	4.5	IO 5A-B	PMA-ZSV45-1	MA-ZSV45-1									
	5.7	IO 5A-C	PMA-ZSV57-1	MA-ZSV57-1									

Price		\$74	\$37	\$37	\$32	\$74	\$74	\$74	\$74	\$31		\$37	
--------------	--	-------------	-------------	-------------	-------------	-------------	-------------	-------------	-------------	-------------	--	-------------	--

Instruments for Scan Body

	Code	Price
Elos Accurate® Scan Body Driver 16mm	C18512	\$73
Elos Accurate® Scan Body Driver 25mm	C13485	\$73
Elos Accurate® Scan Body Driver 92mm	C13048	\$73

Instruments for Print Model Analogs

	Code	Price
Elos Accurate® Analog Pliers	PMA-AIT-1	\$147
Elos Accurate® Analog Insertion Pin	PMA-AIP-2	\$42
Elos Accurate® Analog Insertion Screw	PMA-AIS-5	\$37

Hexlobular Screwdriver

	Code	Price
Elos Accurate® Prosthetic Screwdriver 18mm	PS-AH18-1	\$63
Elos Accurate® Prosthetic Screwdriver 26mm	PS-AH26-1	\$63
Elos Accurate® Prosthetic Screwdriver 34mm	PS-AH34-1	\$63

Torque Wrench and Drivers

	Code	Price
Elos Torque Wrench Kit	C12374	\$315
Elos Driver 4mm Square	C8521	\$79
Elos Driver ISO-1797	C8381	\$79

