

CAMLOG CAD-CAM 3Shape Library Rev F Release Notes

What's New

- Added the new Multi-Use scan bodies for Camlog, Conelog, iSy, and COMFOUR
- Renamed existing Camlog, Conelog, iSy, and COMFOUR scan body files to include "PEEK"

Supported Components



Abutments

- Custom Ti for Camlog, Conelog, iSy
- Ti Base for Camlog, Conelog, iSy, COMFOUR
- Ti Base Bridge for Camlog, Conelog
- Abutment Base for iSy

Scan Bodies

- Authentic: Titanium Scan Bodies, PEEK Scan Bodies, Multifunctional Cap
- 3rd party (Camlog only): Elos, 3Shape, NT Trading

Analogs

- Camlog/Conelog authentic DIM analogs (K/C/J3025.xxxx)

Workflows

- Up to 25° Angled Screw Channel with Ti Bases
- Custom Ti, Temporary & Zirconia restorations

