

BioHorizons CAD-CAM Bars 3Shape Library Rev E Release Notes

What's New

- Updated model retention files – now includes vertical verifier window to help ensure analogs are fully seated in printed model

Supported Components



Connections

- Implant Level: BioH 3.0mm, 3.5mm, 4.5mm, 5.7mm
- Abutment Level: BioHorizons Multi-Unit

○ Implant level:



BioHorizons (PGTSB, PBT SB, PYTSB, TP3TSB)



Elos IO (IO-11A-x)



3Shape (BI3, BIY, BIG, BIB)

Scan Bodies

○ Abutment level:



BioHorizons (PXMUTSB)



Elos IO (IO-11B-A)



Open Implants (3915)



Imetric MU

