

BioHorizons Trimmable Prosthetics 3Shape Library Rev 05 Release Notes

What's New

- Updated model retention files – now includes vertical verifier window to help ensure analogs are fully seated in printed model

Supported Components



Abutments

- Easy Ti Temporary Abutment (Hexed, Non-Hexed)
- PEEK Temporary Cylinder (Hexed, Non-Hexed)
- BioHorizons Multi-Unit Ti Coping

Scan Bodies

- BioHorizons: Titanium Scan bodies*, PEEK Scan Abut, Snap Scan (3 heights), Multi-unit scan caps* (e.g., PXMUCCC)
- 3rd party: Elos*, 3Shape, NT Trading

Analogs

- BioHorizons authentic lab analogs (e.g., PYIA)

Workflows

- Full-Arch Restorations
- Denture Conversion

*NOTE: BioHorizons Multi-Unit Ti Coping only supports indicated scan bodies.

